TMX®
CONNECT

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
TMX®
CONNECT

FREE WEBAPP

- The QR code on the sticky pad you will take you to our free webapp
- Free acute pain exercises
- Free 30-day-pain-free programmes with possibility of individual documentation of progress
- TMX supports you with questions: info@tmx-trigger.de

TMX® Sticky-Pad
INNOVATIVE GEL PAD • STICKS TO ALMOST ANY WALL
REMOVABLE • REUSABLE

HOW TO APPLY THE
ADHESIVE PADS

Underneath our triggers are now the new, innovative **TMX** Sticky-Pads. These allow you to attach the triggers to almost any wall surface.

NOT TO BE USED ON: Wallpaper, painted or plastered walls, sensitively repainted/painted surfaces or similar.

CLEANING

If the adhesive strength of the Sticky Pad wears off, it should be carefully cleaned under lukewarm water (without soap) and then dried in the air.


HOW TO USE
TMX® TRIGGER?

1. Locate your pain area or your trigger point.
2. Place the **TMX** Trigger on the floor / wall and apply proper pressure to the respective area.
3. Take several deep breaths and try to relax despite the pain. Wait maximum of a minute until the initial pain subsides.
IMPORTANT: Make sure to let your muscles loose despite the pressure to ensure a circulation.
4. Now start mobilizing and keep switching back to the holding initial position. Depending on the pain point, you should do this for about 3-5 minutes.
5. After triggering, stretch the corresponding part of the body. You can always add a stretch in between the mobilization and the trigger part.

3 STEPS TO FREE MUSCLES

1. TRIGGER


2. MOBILISE

3. STRETCH

EXAMPLE USES
AT A GLANCE

REACTIONS

Reactions following your self-treatment are normal. Muscle ache, redness on the skin and pain in particular areas after triggering are possible tissue reactions. After using the **TMX** Trigger, you may feel initial worsening after 10-48 hours. Please understand that this is merely a reaction!

All uses must be carried out exactly as shown in the usage images (see reverse). The **TMX** Trigger must be positioned stably on the ground or wall. Always remember your breathing. Do not hold your breath. Keep in mind that the first 2-3 minutes may be painful. Your muscle adapts to this stimulus and the pain is alleviated. Once you have mastered an exercise, you can intensify it. If the intensity and duration of your pain increase considerably, stop the exercise immediately.

INDICATIONS

Muscular tension, Myofascial trigger points, Movement restrictions, Hard tension, Heel spur, Plantar fasciitis, Calf/knee pain, Piriformis syndrome, LWS / - BWS syndrome, Shoulder/neck tension, Impingement syndrome (For further ID see our website)

NOTES

Due to the targeted pressure, the trigger method is painful. After around two minutes, the tissue adjusts to the pressure and relaxes. Please be aware that, despite the pressure, you must completely relax your muscle to ensure tissue circulation.

- The **TMX** Trigger is not a replacement for any physiotherapist, osteopath or doctor
- Not suitable for children aged under 12
- Use of the **TMX** Trigger is at your own risk
- Only use in line with all indications and all contraindications
- Contact a doctor if, after use, pain arises that has not improved after seven days at the latest

CONTRAINDICATIONS

Acute inflammation, redness, swelling, heat, taking anticoagulant medication, acute heart patients, epileptics, fractures, pregnancy (first 3 months and in the last 3 months), thrombosis, severe vascular disease (You can find more AI on our website)

NOTE ON USE: All serious incidents arising in connection with the product are to be reported to the manufacturer and the responsible authorities of the member state in which the user is resident.

GET MORE INFORMATION
WWW.TMX-TRIGGER.DE

THANK YOU FOR YOUR TRUST. DISCOVER THE WHOLE
WORLD OF TMX REGENERATION HERE:

WWW.TMX-CONNECT.DE